

ABSHealthcareLite

QC Engineer Quick Start Guide

MOD-01 — Setup, Commands, and Manual QC Workflow

Version 1.0

Prepared: 19 July 2026

Environment: QC Docker (non-production)

Confidential — ABSHealthcareLite - QC Use Only

Table of Contents

1. Documents you need
2. QC environment URLs
3. QC server — deploy and verify
4. Local development (optional)
5. Manual QC workflow
6. Quick command reference
7. Roles during testing
8. Important rules
9. Support references
10. Contact and escalation

QC Engineer Quick Start Guide

ABSHealthcareLite Pilot — MOD-01

Document control	Value
---	---
Document version	1.0
Prepared date	19 July 2026
Environment	QC Docker (non-production)
Intended audience	Independent QC engineers
Classification	Internal QC working document

****Credential handling:**** Passwords are ****not printed in this document****. Obtain credentials from the secure channel or ``local/001-QC-Credentials.txt`` on the QC workstation (never email or chat).

1. Documents you need

Document	Path / location	Purpose
---	---	---
Manual QC test guide (PDF)	<code>`docs/AI-QC/manual-qc/pdf/001-CompanyTenant-Manual-QC-v1.0.pdf`</code>	Full test cases (142 cases)
Result template	<code>`docs/AI-QC/manual-qc/results/001-CompanyTenant-Manual-QC-Result-Template.md`</code>	Copy and fill during execution
QC credentials	<code>`local/001-QC-Credentials.txt`</code> (secure channel only)	Host and tenant login passwords
Docker QC guide	<code>`docs/QC/ABSHealthcareLite-Pilot-Docker-QC-Guide.md`</code>	Server and container details
Deployment report	<code>`docs/AI-QC/manual-qc/evidence/001-CompanyTenant/deployment/`</code>	Latest deploy status

****Bangla (বাংলা):**** প্রথমে Manual QC PDF খুলুন, result template কপি করুন, credentials secure channel থেকে নিন। Test case গুলো PDF অনুযায়ী চালান।

2. QC environment URLs

Item	Value
---	---
Application	<code>`http://192.168.2.44:3000`</code>
Health check	<code>`GET http://192.168.2.44:3000/api/health`</code>
Host login	<code>`http://192.168.2.44:3000/host/login`</code>
Tenant login	<code>`http://192.168.2.44:3000/login`</code>
Host username	<code>`admin.abs`</code>
Tenant username	<code>`laila.hasan`</code>
Passwords	Provided separately through secure channel
App container	<code>`abs-healthcare-pilot-app`</code>
Database container	<code>`abs-healthcare-pilot-db`</code>
Database name	<code>`abs_healthcare_pilot`</code>

Expected health response: JSON with `"status": "ok"` and `"database": "connected"`.

3. QC server — deploy and verify (operator)

Run on the QC Docker host after ``git pull``:

```
cd abs-healthcare-pilot
git pull origin main
docker compose up -d --build
docker compose exec app npx prisma db seed
```

Verify:

```
curl http://192.168.2.44:3000/api/health
docker ps --filter name=abs-healthcare-pilot
```

Note: Container startup runs migrations only. Seed is **manual** (command above).

4. Local development (optional — for developers)

Use only if testing on a developer workstation, not for official QC sign-off.

```
cd abs-healthcare-pilot
copy .env.example .env
npm install
npm run db:migrate
npm run db:seed
npm run dev
```

Local URL	Path
---	---
App	<code>`http://localhost:3000`</code>
Health	<code>`http://localhost:3000/api/health`</code>
Host login	<code>`http://localhost:3000/host/login`</code>
Tenant login	<code>`http://localhost:3000/login`</code>

Optional verification:

```
npx tsx scripts/verify-mod01-auth.ts
npm run build
```

5. Manual QC workflow (QC engineer)

Step 1 — Pre-test checks

1. Open health endpoint and save the JSON response as evidence.
2. Confirm Docker containers are running.
3. Use a clean or isolated browser profile.
4. Record tester name, date, and environment in the result file.
5. Do **not** modify other testers' data in shared QC.

Step 2 — Prepare result file

Copy the result template to a dated file, for example:

``001-CompanyTenant-Manual-QC-Result-20260719.md``

Step 3 — Execute test cases

1. Open ``001-CompanyTenant-Manual-QC-v1.0.pdf``.

- 2. Execute every `MOD01-\*` test case in order (unless a dependency note says otherwise).
- 3. For each case record: **Actual Result**, **Status**, **Evidence**, **Tester Remarks**.
- 4. Status values: `PASS`, `FAIL`, `BLOCKED`, `NOT APPLICABLE`, or `NOT RUN`.

Step 4 — Test data naming

Data item	Convention
---	---
Tenant code	`QC01-YYYYMMDD-HHMM`
Tenant name	`MOD01 QC Tenant`
Branch code	`QC-BR-01`
Unique suffix	Execution timestamp, e.g. `20260719-1958`

Use a **new timestamp suffix** for each run. Log every tenant ID and code created.

Step 5 — Save evidence

Evidence type	Folder
---	---
Screenshots	`docs/AI-QC/manual-qc/evidence/001-CompanyTenant/screenshots/`
Database output	`docs/AI-QC/manual-qc/evidence/001-CompanyTenant/database/`
Security tests	`docs/AI-QC/manual-qc/evidence/001-CompanyTenant/security/`

Screenshot naming: `MOD01-<TCID>-StepNN.png` (example: `MOD01-AUTH-001-Step01.png`).

Step 6 — Sign-off

Complete sign-off only when exit criteria in Manual QC guide Section 13 are met. Do not mark UAT or production approved from this pack alone.

6. Quick command reference

Task	Command
---	---
Pull latest code	`git pull origin main`
QC Docker deploy	`docker compose up -d --build`
QC database seed	`docker compose exec app npx prisma db seed`
Health check	`curl http://192.168.2.44:3000/api/health`
List containers	`docker ps --filter name=abs-healthcare-pilot`
MOD-01 auth verify (local)	`npx tsx scripts/verify-mod01-auth.ts`
Regenerate Manual QC PDF	`node scripts/generate-manual-qc-pdf.mjs`
Local dev server	`npm run dev`

7. Roles during testing

Role	Purpose
---	---
Host Admin	Create, edit, activate, suspend, and review tenants
Tenant Admin	Manage assigned tenant only
Unauthorized / Anonymous	Validate login and access rejection
Database Observer	Run approved **read-only** SQL only

8. Important rules

1. **QC only** — target is `192.168.2.44` / `abs\_healthcare\_pilot`, not live production.
2. **No passwords in evidence** — blur or omit password fields in screenshots.
3. **No destructive SQL** — read-only database checks only.
4. **Stop on FAIL** — log defect ID before dependent cases if a blocker is found.
5. **Independent execution** — AI Re-QC PASS does not replace manual test execution.

9. Support references

Topic	Document
---	---
Full test cases	`001-CompanyTenant-Manual-QC-v1.0.pdf`
AI Re-QC baseline	`docs/AI-QC/reports/001-Foundation-ReQC-01.md`
Module specification	`docs/01-CompanyTenant/CompanyTenant_v1.md`
Git repository	`https://github.com/abidaiub/abs-healthcare-pilot.git`
Latest commit	`main` branch — check `git log -1` on QC host before testing

10. Contact and escalation

Situation	Action
---	---
Health check fails	Escalate to deployment operator; do not start functional QC
Credentials missing	Request via secure channel; do not use placeholder passwords
Environment blocked	Mark affected cases `BLOCKED` and reference deployment report
Defect found	Record in result file; attach screenshot evidence